




### PROFILE



For liners with a progressive profile, the size of the accordion distal umbrella changes with the size of the liner.

### SIZING GUIDE

Please refer to the illustrations for recommendations for fitting Alpha Silicone Liners.




### TIP

If the limb circumferences fall into the range of both the Medium Plus 2 and the Large 1:

- choose **Medium Plus 2** if the distal end is bony and tapered
- choose **Large 1** if the distal end is broad and cylindrical

### RETROFITTING

Patients wearing Alpha Hybrid Liners with a progressive profile can be retrofit with Alpha Silicone Liners with a progressive profile in the same base size (Medium Plus to Medium Plus, or Large to Large).



\*A Medium Plus 3 size is available upon request to retrofit patients that are near the upper limit of the Medium Plus size range.

### ITEM NUMBERS

#### Locking

Item #	Size
S352-6590-1	Medium 1
S352-6590-2	Medium 2
S352-6594-1	Medium Plus 1
S352-6594-2	Medium Plus 2
S352-6594-3	Medium Plus 3
S352-6596-1	Large 1
S352-6596-2	Large 2

#### Cushion

Item #	Size
S350-6590-1	Medium 1
S350-6590-2	Medium 2
S350-6594-1	Medium Plus 1
S350-6594-2	Medium Plus 2
S350-6594-3	Medium Plus 3
S350-6596-1	Large 1
S350-6596-2	Large 2

