

Alpha® Silicone Liners

WillowWood is excited to announce its new Alpha Silicone Liner! Our Silicone Liner uses a proprietary, platinum-cured, medical grade silicone designed for comfort and performance. Our unique blend of silicone, Vitamin E, and skin conditioners results in a pleasing surface that is non-greasy, non-tacky, and comfortable against the skin. The Alpha Silicone Liner uses Select fabric that controls pistoning without using a distal matrix. The Select fabric allows for circumferential stretch making it easy to don for numerous limb shapes. The flexible knee panel provides enhanced knee flexion and a broad range of knee motion. The Alpha Silicone Liner may be retrofit with Alpha Hybrid Liners with a progressive profile.

Release Date:
 November 1, 2011

Product Highlights:

- Delivers a Comfortable Fit Using a Medical-Grade Silicone.
- Silicone Formulated to have a Non-Greasy and Non-Tacky Surface.
- Controls Pistoning without a Distal Matrix.
- Flexible Knee Panel Provides Enhanced Knee Flexion
- Retrofits with Thermoplastic Elastomer Liners with Progressive Profile
- Available in Locking and Cusion
- Features a 12-month warranty



Styles: Progressive
 Gel: Silicone
 Colors: Taupe/Gray
 Thickness: 9 mm distal, 6 mm tibial crest, 2.5 mm proximally (Progressive)
 Size: Medium Plus 1, Medium Plus 2, Large 1, Large 2

